User-friendliness and ease of use

The **SUPPORTAIR®** combines the user-friendliness and ease of use of the AIROX HomeCare ventilators with the technical requirements that are made compulsory with taking care of patients in instable respiratory phase at hospital.

The **preference** and the **curve setup menus** that group together on a unique and same display all the adjustable components allow the user to set up its ventilator following its habits.



PREFERENCES

D Storm Serv. YES

Communication software

The **AIROX COMMUNICATION** software is compatible with the **SUPPORTAIR**® but also with the others AIROX HomeCare ventilators.

Various thumbs are available allowing retrieval, edition and print:

- Of patient compliance
- Of events memory summarizing data stocked in the machine such as ventilation alarms, technical alarms or also events linked to the device use or settings.
- Of monitoring providing curves, loops, monitored values and alarms either in real-time or for 24 hours.
 - Of maintenance menu allowing the ventilator software update.

The **SUPPORTAIR®**

is a hospital ventilator that benefits from AIROX HomeCare ventilation know-how regarding the ventilation comfort, ergonomics, reduced maintenance

> but overall in terms of user-friendliness and reliability



Maintenance

The **SUPPORTAIR**[®] benefits from the AIROX know-how in HomeCare ventilation, area that imposes development of ventilators with very high reliability, requiring minimum preventive maintenance and very short intervention times.

Soft: LH010000

p Sensors Calibration

o Insp Flow

o Exhal Flow

o 02 Pressure

D Turbine Speed : -

e 02 Flow

o Patient Pressure

e Valve Pressure

■ Defaults

□ Measures Check

It is characterized by:

- Quick access to the different technical components
- Very few components
- Micro-turbine with no maintenance required for a period of 15.000 hours

- RS 232 port allowing event retrieval, monitoring as well as new software's download.
- Maintenance menu allowing the following:

0.0mbar / 0.0mbar

0.0mbar

134.0 lpm

13.5 lpm

0.0 lpm

0.0 bar

- MAINTENANCE as software version, machine hours meter, anomalies
 - Measurements and direct controls (voltage, temperature,...) allowing making quick diagnostics,
 - 3 Sensors calibration (flow, pressure, FiO₂),

Technical features

O₂ connexion

Class

 Dimensions (H,L,P)
 H = 154 mm, W = 235 mm, D = 315 mm

 Weight
 4,9 kg

 Power supply
 115 / 230 V ± 10% - 50 / 60 Hz

 Electric consumption
 90 VA max

Battery (type, autonomy) 25,2 V = 4,4 Ah Lithium-Ion type with quick charge, autonomy up to 11 hours depending on settings.

High pressure O₂ inlet: 280-600 kPa, (2,8 – 6 bars)

Device and accessoires references

Reference	List of items				
4096200 (Global) 4096200 (UK)	$\label{eq:SUPPORTAIR} Supplied with 230 V power supply cable, O_2 supply tubing, FiO_2 measurement kit, user manual, single use patient circuit with valve, SpO_2 \text{ sensor with its extension.}$				
Accessories					
2739200 5092300 3816000 4096000 3816500 4096400 4096500 4095600 3809999 3823099 3814100	O ₂ supply tubing (3m length) Reusable sensor for pulse oxymeter (SpO ₂), Extension cable for SpO ₂ sensor Remote alarm Remote alarm cable Trolley Patient circuits holder AIROX COMMUNICATION software pack Reusable exhalation block Single patient use exhalation stock FiO ₂ measurement kit (T adaptor, FiO ₂ cell, connexion cable, deflector)				
Disposable					
5092100 5092200 2963300 2964200	Single use, single branch patient circuit with exhalation valve Single use, double branch patient circuit with exhalation valve Air inlet combi filter, bag of 6 (foam + thin particle) FiO ₂ Cell				





PARC D'ACTIVITÉS PAU-PYRÉNÉES - L'ÉCHANGEUR B.P. 833 - 64 008 PAU cedex - FRANCE

Tel: +33 (0) 5 59 14 02 02 Web Site: www.airox.fr Fax: +33 (0) 5 59 14 02 00 E-mail: contact@airox.fr SUPPORTAIR®

Titration platform for HomeCare ventilation

Titration platform

The SUPPORTAIR® is a titration platform intended for hospitalized patients suffering from Chronic Respiratory Failures (CRF). Its polyvalence and numerous functions allow it to take care of those chronic patients as soon as they are admitted at hospital for acute exacerbation, all along their health state evolution, and until they are stabilized and ready to go back home.

Thus, the SUPPORTAIR® finds its place in all units like intensive care (ICU), respiratory assistance units. continuous care and observation units, pneumology department,...

Settings portability

Settings and parameters defined and validated step by step with the SUPPORTAIR® (ventilation titration) are directly "portable" onto the whole AIROX HomeCare ventilators range thanks to a pneumatic system and ergonomics voluntary common to all.

This settings portability feature, unique and specific to the SUPPORTAIR®, considerably makes easier the prescription phase.



Polyvalency, exhaustive ventilatory functions and monitoring

The SUPPORTAIR® is a mixed ventilator with large poly-

Ilt offers the following functionalities:


- Non invasive and invasive ventilation
- Leakage and valve ventilation
- Volumetric (CV/ACV/SIMV) and pressure support (CPAP/PSV-S/PSV-ST/PCV/PACV) ventilation
- Adults and children (>5kg) ventilation
- Ventilation with single or double branch circuits
- Oxygen blender allowing FiO₂ adjustment from 21 to
- Complete monitoring (pressure/flow curves, loops, SpO_{2....})
- Exhaustive ventilatory functions (Vt TARGET, Ti min/max, 100% O_{2,...})
- Internal battery providing up to 11 hours autonomy

O₂ integrated blender

The SUPPORTAIR® is equipped with a high pressure oxygen inlet (wall oxygen) and an integrated blender that allows adjusting FiO2 from 21 to 100%. A specific key (25%) allows ram O2 induction at 100% for 2 minutes (suction,...).

Moreover, FiO₂ monitoring and associated alarms configurations are available directly on the ventilator





Monitoring display

either via curves, loops or nu-meric values. This display allows getting in glance a complete view o the on going ventilation.





SpO₂ integrated monitoring

AIROX

During ventilation titration, it is necessary to watch out oxygenation in

available in the SUPPORTAIR®

ventilation parameters.

a non invasive way. The FiO2 setting and

monitoring association together with the ventilator SpO₂ integrated surveillance allows optimization of

the patients oxygenation management. The SpO₂

makes possible instantaneous or recorded analysis

and review of trends correlated with the others

ventilation functions gear is loaded on the SUPPORTAIR® in order to make easier the patients management particularly during ventilation instability periods as well as during the titration phase

A complete

The Vt TARGET function enables maintaining target volume (Vt) chosen by the clinician, thanks to pressure automatic adjustment into a predetermined margin even in case of nocturnal compliance variation (change of position or of sleep phase....), or of resistance variation

A better Vt and minute volume stability all along the night generating a better comfort, a better alveolar ventilation optimization and a potential reduction of the nocturnal hypoventilations are part of advantages provided by that function.


Vt TARGET

(secretions, aperture of the upper airways,...)

Adjustable ventilation parameters

The SUPPORTAIR® offers a parameters setting menu specific to each mode whether it is active or not, composed of list of adjustable parameters, a pressure bar graph and a monitoring window in order to set up the ventilator in the best conditions.






Summary of adjustable ventilation parameters


Ventilation modes	CPAP	PSV-S / PSV-ST	PCV / PACV	CV/ACV	SIMV
Tidal volume (Vt)	X	X	X	50 - 2000 ml	50 - 2000 m
Vt TARGET	X	50 - 2000 ml	50 -2000 ml	X	X
Maxi P	X	8 - 60 mbar	8 - 60 mbar	X	X
Sigh	X	X	X	yes/no	X
Sigh Vt	X	X	X	(1.0 - 2.0) x Vt	X
Sigh rate	X	X	X	50-250 cycles	X
Inspiratory pressure (IPAP)	X	5 - 60 mbar	5 - 60 mbar	X	5 - 40 mbar
Expiratory pressure (EPAP)	X	0 - 20 mbar	0 - 20 mbar	0 - 20 mbar	0 - 20 mbar
CPAP	4 - 20 mbar	X	X	X	X
Rise Ti	X	1 - 4	1 - 4	X	X
Rate	X	4 - 40 bpm	5 - 60 bpm	5 – 60 bpm	5 - 40 bpm
I/E ratio	X	X	1.0 - 3.0	1.0 - 3.0	1.0 - 3.0
I/T ratio	X	X	25 - 50 %	25 - 50 %	25 - 50 %
Ti min	X	0.1 - 2,8 s	X	X	X
Ti max	X	0.8 - 3 s	X	X	X
Trigger I	X	1(Pedia) - 5	1(Pedia) - 5	1(Pedia) - 5	1(Pedia) - 5
Trigger E	X	(-5)-(-95)%	X	X	X
FiO ₂	21 - 100 %	21 - 100 %	21 - 100 %	21 - 100 %	21 - 100 %
Apnea	X	3-30 s	X	X	3-30 s
Ramp	X	X	X	R/D	X
SIMV rate	X	Х	Х	Х	2 - 20 bpm

Adjustable alarms and monitoring

Each mode presents a specific monitoring and alarm menu providing:

Display of patient values





- monitored during the cycle Setting of high and low alarm parameters
- Additional window of monitoring and alarm message(s) display
- Alarms memory